

Fleet POCT Tube Cheat Sheet

🩸 All tests use whole blood

Test / Cartridge	Device	Tube	Additive	Notes
CHEM8+	i-STAT	Green Top	Lithium Heparin	Test within 10 mins of collection.
CG4+	i-STAT	Green Top	Lithium Heparin	Test immediately.
cTnl (Troponin I)	i-STAT	Green Top	Lithium Heparin	Test within 30 mins of collection.
Troponin I (POCT)	POCT	Green Top	Lithium Heparin	3 drops of whole blood. Read at 15 mins.
PT/INR	i-STAT	Red Top	No additive	Fingerstick preferred. Test immediately.
Mono Test	POCT	Lavender	EDTA	2 drops of whole blood and then 1 drop of R1 buffer. Read at 5 mins.



BD 367880



BD 367861



BD 366408