

**Appendix 1 to
Annex D to
CF H Svcs Gp Instruction 4080-02**

Patients and QC log sheet - Rapid Strep group A Test Kit

Materials Required:

1. QuickVue In-Line Strep A Test, NSN: 6550-21-870-8620
2. Positive/Negative Control Set, NSN: 6550-CF-002-9346
3. Timer

Method Outline: Refer to SOP & package insert for testing procedures

1. Perform Quality Control each time a new box is open and thereafter once a month using external QC swabs provided.
2. Read and record all results on log sheet and on the patient's requisition.
3. Monitor Built in Controls and record results on log sheet.
4. Lab recommends confirmatory test (throat swab) send to the lab.
5. Submit log sheets to laboratory's Sgt or to the regional WO at the end of each quarter.

Month:

Year:

Unit:

Reviewed by:

Kit lot #:		Open date:		Exp. date:		(do not use after exp. date)			
POS QC lot # :		Open date:		Exp. date:		(do not use after exp. date)			
NEG QC lot #:		Open date:		Exp. date:		(do not use after exp. date)			
Date	Last Name	First Name	Svc #	Result pos/neg	Performed by	Scan into CFHIS	Built-in Controls		
							Green color Yes/No	Blue line QC Yes/No	Result (Test line) White or Pink
	Positive Control Swab								
	Negative Control Swab								

Month: **Year:** **Unit:** **Reviewed by:**

PROTECTED B (When completed)